

Implant Component Summary 2024.V3

December 12, 2024

Brand	System	Level	Platform	Internal code	Abutment ¹ w/ Ti-Base	Bridge ¹ w/ Ti-Base	Torque Ncm	Angle	Digitek ² Abutment	Digitek ² Bar/ Bridge	Original Custom Milling Abutment	Torque Ncm	Angle	MTX ³ Zir	Screw torque Ncm	MTX Ti / CoCr ²		Torque Ncm	Angle	PEEK ⁵ Customized Healing	
																Osteon	Digitek				
Nobel Biocare	Replace (NR)	Fixture	NP	NR-100	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	25	40°	TBC	
			RP	NR-101	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	25	40°	TBC	
			WP	NR-102	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	25	40°	TBC	
	Active/Conical (NA)	Fixture	6.0	NR-103	✓	✓	35	25°	✓	✓	—	35	40°	✓	—	✓	✓	25	40°	TBC	
			3.0	NA-100	✓	✓	15	20°	✓	✓	—	15	30°	—	—	—	—	—	—	TBC	
			NP	NA-101	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	25	40°	TBC	
	Branemark (NB)	Fixture	RP	NA-102	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	25	40°	TBC	
			WP	NA-103	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	25	40°	TBC	
			NP	NB-100	✓	✓	35	20°	✓	✓	—	35	30°	✓	20	✓	✓	20	30°	TBC	
	MUA (NU)	Abutment	RP	NB-101	✓	✓	35	25°	✓	✓	—	35	40°	✓	25	✓	✓	30	40°	—	
			WP	NB-102	✓	✓	35	25°	✓	✓	—	35	40°	✓	—	✓	✓	25	40°	—	
	On1 (NO)	Abutment	RP	NU-100	✓	✓	15	40°	✓	✓	—	15	40°	✓	—	✓	✓	15	40°	—	
			WP	NU-101	✓	✓	15	40°	✓	✓	—	15	40°	✓	—	✓	✓	25	40°	—	
	Straumann	Tissue Level (ST)	Fixture	NNC	ST-100	✓	✓	35	20°	✓	✓	—	35	30°	✓	—	—	—	—	—	TBC
RN				ST-101	✓	✓	35	20°	✓	✓	✓*	35	40°	✓	25	✓	✓	25	40°	TBC	
WN				ST-102	✓	✓	35	25°	✓	✓	✓*	35	40°	✓	—	✓	✓	25	40°	TBC	
Bone Level (SB)		Fixture	SC	SB-109	—	—	—	—	✓	✓	—	15	—	—	—	—	—	—	—	TBC	
			NC	SB-100	✓	✓	35	20°	✓	✓	✓*	35	30°	✓	20	✓	✓	25	30°	TBC	
			RC	SB-101	✓	✓	35	20°	✓	✓	✓*	35	30°	✓	20	✓	✓	25	30°	TBC	
Screw-Retained Abutments (SBU) BLX (SX)		Abutment	3.5	SBU-100	—	—	—	—	✓	✓	—	15	40°	—	—	✓	✓	15	40°	—	
			4.6	SBU-101	—	—	—	—	✓	✓	—	15	40°	—	—	✓	✓	15	40°	—	
Tissue Level X (TLX)		Fixture	RB/ WB	SX-102	—	—	—	—	✓	✓	✓*	35	30°	—	—	—	✓*	35	30°	TBC	
			WB	SX-103	—	—	—	—	✓	✓	✓*	35	30°	—	—	—	✓*	35	30°	TBC	
			NT	TLX-100	✓	✓	35	35°	✓	✓	✓*	35	35°	—	—	—	✓	35	35°	TBC	
Megagen		Anvridge (MA)	Fixture	—	MA-100	✓	✓	35	10°	✓	✓	—	35	10°	—	—	—	—	—	—	TBC
				3.0	DA-100	✓	✓	15	15°	✓	✓	—	15	20°	—	—	—	—	—	—	TBC
				3.5/4.0	DA-101	✓	✓	25	15°	✓	✓	—	20	20°	✓	20	✓	✓	25	20°	TBC
	UniAbutment (DAU)	Abutment	4.5/5.0	DA-102	✓	✓	25	15°	✓	✓	—	25	20°	✓	25	✓	✓	25	20°	TBC	
			2.0°	DAU-100	—	—	—	—	✓	✓	—	15	40°	—	—	✓	✓	15	40°	—	
			4.5°	DAU-101	—	—	—	—	✓	✓	—	15	60°	—	—	✓	✓	15	60°	—	
	AstraTech EV (DV)	Fixture	3.6	DV-101	✓	✓	25	20°	✓	✓	—	25	20°	✓	20	✓	✓	25	20°	TBC	
			4.2	DV-102	✓	✓	25	20°	✓	✓	—	25	20°	✓	25	✓	✓	25	20°	TBC	
			4.8	DV-103	✓	✓	25	20°	✓	✓	—	25	20°	✓	—	✓	✓	25	20°	TBC	
	UniAbutment EV (DVU) Multi base(DVM) FrialitXIVE (DX)	Abutment	5.4	DV-104	—	—	—	—	✓	✓	—	25	15°	—	—	✓	✓	25	15°	TBC	
			Straight	DVU-100	—	✓	15	40°	—	✓	—	15	40°	—	—	✓	✓	25	40°	—	
	Ankylos (DY)	Fixture	3.4	DX-100	✓	✓	25	25°	✓	✓	—	24	40°	✓	20	✓	✓	25	40°	TBC	
			3.8	DX-101	✓	✓	25	25°	✓	✓	—	24	40°	✓	20	✓	✓	25	40°	TBC	
			4.5	DX-102	✓	✓	25	25°	✓	✓	—	24	40°	✓	—	✓	✓	25	40°	TBC	
5.5			DX-103	✓	✓	25	25°	✓	✓	—	24	40°	✓	—	✓	✓	25	40°	TBC		
C/X			DY-100	✓	✓	15/25	15°	—	✓	—	15	10°	—	—	—	—	—	—	—	TBC	
Zimmer	Tapered Screw-Vent (ZT)	Fixture	3.5	ZT-100	✓	✓	30	60°	✓	✓	—	30	60°	✓	25	✓	✓	25	60°	TBC	
			4.5	ZT-101	✓	✓	30	60°	✓	✓	—	30	60°	✓	25	✓	✓	25	60°	TBC	
			5.7	ZT-102	✓	✓	30	60°	✓	✓	—	30	60°	✓	—	✓	✓	25	60°	TBC	
3i Biomet	Certain (BC)	Fixture	3.4	BC-100	✓	✓	20	25°	✓	✓	—	20	40°	✓	25	✓	✓	25	40°	TBC	
			4.1	BC-101	✓	✓	20	25°	✓	✓	—	20	40°	✓	25	✓	✓	25	40°	TBC	
			5.0	BC-102	✓	✓	20	25°	✓	✓	—	20	40°	✓	—	✓	✓	25	40°	TBC	
External Hex (BO)	Fixture	3.4	BO-100	✓	—	35	20°	✓	✓	—	20	30°	—	—	✓	✓	30	30°	—		
		4.1	BO-101	✓	—	35	25°	✓	✓	—	20	40°	—	—	✓	✓	30	40°	—		
		5.0	BO-102	✓	—	35	25°	✓	✓	—	20	40°	—	—	✓	✓	30	40°	—		
Osstem	TS (OT)	Fixture	Mini	OT-100	✓	✓	20	20°	✓	✓	—	20	20°	✓	20	✓	✓	25	20°	TBC	
			Standard	OT-101	✓	✓	30	20°	✓	✓	—	30	20°	✓	25	✓	✓	25	20°	TBC	
Camlog	Camlog (CC)	Fixture	3.3	CC-100	✓	✓	20	25°	✓*	✓	✓*	20	40°	—	—	—	—	—	—	TBC	
			3.8	CC-101	✓	✓	20	25°	✓	✓	✓*	20	40°	—	—	—	—	—	—	TBC	
			4.3	CC-102	✓	✓	20	25°	✓	✓	✓*	20	40°	—	—	—	—	—	—	TBC	
Conelog (CO)	Fixture	5.0	CC-103	✓	—	20	25°	✓	✓	✓*	20	40°	—	—	—	—	—	—	—	TBC	
		3.3	CO-100	✓	✓	20	25°	✓*	✓	✓*	20	40°	—	—	—	—	—	—	TBC		
		3.8	CO-101	✓	✓	20	25°	✓	✓	✓*	20	40°	—	—	—	—	—	—	TBC		
ISY (CI)	Fixture	4.3	CO-102	✓	✓	20	25°	✓	✓	✓*	20	40°	—	—	—	—	—	—	—	TBC	
		5.0	CO-103	✓	✓	20	25°	✓	✓	✓*	20	40°	—	—	—	—	—	—	—	TBC	
		—	CI-100	✓	✓	—	—	—	—	✓	✓*	20	40°	—	—	—	—	—	—	—	TBC
Anthogyr	Axiom BL (AA)	Fixture	—	AA-100	✓	✓	25	15°	✓	✓	—	25	10°	—	—	✓	✓	20	10°	TBC	
			4.0	AT-100	—	—	—	—	✓	✓	—	25	40°	—	—	—	—	—	—	—	
	Multi Unit (AAU)	Abutment	4.8	AT-101	—	—	—	—	✓	✓	—	25	40°	—	—	—	—	—	—	—	
			Narrow	AAU-100	—	—	—	—	—	✓	—	15	40°	—	—	—	—	—	—	—	
Biotech	Kontakt (TK)	Abutment	Standard	AAU-101	—	—	—	—	—	✓	15	40°	—	—	—	—	—	—	—	—	
			Narrow	AI-100	—	Hold	25	60°	—	Hold	—	25	60°	Hold	25	—	Hold	25	60°	—	
			Regular	AI-101	—	Hold	25	60°	—	Hold	—	25	60°	Hold	25	—	Hold	25	60°	—	
BTI	Interna (TI)	Abutment	TK-100	TKU-100	✓	✓	20	10°	✓	✓	—	20	10°	—	—	—	—	—	—	TBC	
			Standard	TKU-101	—	—	—	—	✓	✓	—	20	40°	—	—	—	—	—	—	—	
			4.1	TI-101	✓	✓	35	40°	✓	✓	—	35	40°	—	—	—	—	—	—	—	
Medentis	ICX (MEI)	Fixture	Straight	TIU-100	—	—	—	—	✓	✓	—	20	40°	—	—	—	—	—	—	—	
			Angulado	TIU-101	—	—	—	—	✓	✓	—	20	40°	—	—	—	—	—	—	—	
			MEI-100	MEI-100	✓	✓	30	20°	✓	✓	—	30	20°	—	—	—	—	—	—	—	
MIS	CI (MIC)	Fixture	NP	MIC-100	✓	✓	30	12°	✓*	✓	—	30	40°	—	—	✓	✓	25	40°	TBC	
			SP	MIC-101	✓	✓	30	12°	✓*	✓	—	30	40°	—	—	✓	✓	25	40°	TBC	
			WP	MIC-102	✓	✓	30														